



3.5mm Diamond Tipped Burr 35 Cms Shaft (Imported German Made)

3.5mm Diamond Tipped Burr 35 Cms Shaft (Imported German Made)

[Read More](#)

SKU: MDB350

Price:

Stock: instock

Categories: [Drills, Perforators, Trephines, Endoscopic Transforaminal, Neuro & Spine, Perforators, Trephines](#)

Tags: [MDB350](#)

Product Description

3.5mm Diamond Tipped Burr 35 Cms Shaft (Imported German Made)