



Curved  
46.02.12 30.4 cm - 12"  
46.02.16 40.6 cm - 16"  
46.02.18 45.7 cm - 18"  
46.02.20 50.8 cm - 20"  
46.02.24 60.9 cm - 24"

# Anderson - Brand Shunt Tunneler Curved 45.7 Cm - 18"

Anderson - Brand Shunt Tunneler Curved 45.7 Cm - 18"

[Read More](#)

**SKU:** OTN030

**Price:**

**Stock:** instock

**Categories:** [Ortho](#), [Tunneler](#)

**Tags:** [46.02.18](#)

## Product Description

Anderson - Brand Shunt Tunneler Curved 45.7 Cm - 18"