



Curved
46.02.12 30.4 cm - 12"
46.02.16 40.6 cm - 16"
46.02.18 45.7 cm - 18"
46.02.20 50.8 cm - 20"
46.02.24 60.9 cm - 24"

Anderson - Brand Shunt Tunneler Curved 45.7 Cm - 18"

Anderson - Brand Shunt Tunneler Curved 45.7 Cm - 18"

[Read More](#)

SKU: OTN030

Price:

Stock: instock

Categories: [Ortho](#), [Tunneler](#)

Tags: [46.02.18](#)

Product Description

Anderson - Brand Shunt Tunneler Curved 45.7 Cm - 18"