



Straight
46.01.12 30.4 cm - 12"
46.01.16 40.6 cm - 16"
46.01.18 45.7 cm - 18"
46.01.20 50.8 cm - 20"
46.01.24 60.9 cm - 24"

Anderson - Brand
Shunt Tunneler

Anderson - Brand Shunt Tunnelers Straight 40.6 Cm - 16"

Anderson - Brand Shunt Tunnelers Straight 40.6 Cm - 16"

[Read More](#)

SKU: OTN010

Price:

Stock: instock

Categories: [Ortho](#), [Tunnelers](#)

Tags: [46.01.16](#)

Product Description

Anderson - Brand Shunt Tunnelers Straight 40.6 Cm - 16"