



**Straight**  
46.01.12 30.4 cm - 12"  
46.01.16 40.6 cm - 16"  
46.01.18 45.7 cm - 18"  
46.01.20 50.8 cm - 20"  
46.01.24 60.9 cm - 24"

Anderson - Brand  
Shunt Tunneler

# Anderson - Brand Shunt Tunneler Straight 45.7 Cm - 18"

Anderson - Brand Shunt Tunneler Straight 45.7 Cm - 18"

[Read More](#)

**SKU:** OTN015

**Price:**

**Stock:** instock

**Categories:** [Ortho](#), [Tunneler](#)

**Tags:** [46.01.18](#)

## Product Description

Anderson - Brand Shunt Tunneler Straight 45.7 Cm - 18"