



# Autoclavable Container With Silicon Inserts & Bridge To Hold 2 Endoscopes

Autoclavable Container With Silicon Inserts & Bridge To Hold 2 Endoscopes

[Read More](#)

**SKU:** DMB403

**Price:**

**Stock:** instock

**Categories:** [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

**Tags:** [DMB403](#), [FL108](#)

## Product Description

Autoclavable Container With Silicon Inserts & Bridge To Hold 2 Endoscopes