



# Autoclavable Container With Silicon Inserts & Bridge To Hold Tfe Set

Autoclavable Container With Silicon Inserts & Bridge To Hold Tfe Set

[Read More](#)

**SKU:** DMB401

**Price:**

**Stock:** instock

**Categories:** [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

**Tags:** [DMB401](#), [FL104](#)

## Product Description

Autoclavable Container With Silicon Inserts & Bridge To Hold Tfe Set

---