



Dilator 6.3mm 24 Cms Id 1.5mm Single Bore Taper Tip 6.3mm Chinese Scope

Dilator 6.3mm 24 Cms Id 1.5mm Single Bore Taper Tip 6.3mm Chinese Scope

[Read More](#)

SKU: DMB405A

Price:

Stock: instock

Categories: [Dilator](#), [Endoscopic Transforaminal](#), [Neuro & Spine](#)

Tags: [DMB405A](#), [FL186](#)

Product Description

Dilator 6.3mm 24 Cms Id 1.5mm Single Bore Taper Tip 6.3mm Chinese Scope