



DissectorCvd Fibre Handle 3mm Tip Working Length 30cms

DissectorCvd Fibre Handle 3mm Tip Working Length 30cms

[Read More](#)

SKU: DMB517

Price:

Stock: instock

Categories: [Dissectors, Raspatory, Elevators etc](#), [Elevators etc](#),
[Endoscopic Interlaminar, Neuro & Spine, Raspatory](#)

Tags: [DMB517](#), [FL526](#)

Product Description

DissectorCvd Fibre Handle 3mm Tip Working Length 30cms