



Endoscopic Tubular Kerrison 3.3mm 90 Deg Up Rotatable

Endoscopic Tubular Kerrison 3.3mm 90 Deg Up Rotatable

[Read More](#)

SKU: DMB210R

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB210R](#)

Product Description

Endoscopic Tubular Kerrison 3.3mm 90 Deg Up Rotatable