



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 2mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 2mm

[Read More](#)

SKU: DMB529

Price:

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB529](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 2mm