



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 3mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 3mm

[Read More](#)

SKU: DMB531

Price:

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB531](#), [FL573](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 3mm
