



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 4mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 4mm

[Read More](#)

SKU: DMB533

Price:

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Rongeurs](#), [Neuro & Spine](#)

Tags: [DMB533](#), [FL575](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 4mm
