



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 5mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 5mm

[Read More](#)

SKU: DMB535

Price:

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB535](#), [FL576](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 5mm
