



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up3mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up3mm

[Read More](#)

SKU: DMB537

Price:

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB537](#), [FL577](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up3mm
