



Needle Holder Mathieu-Kocher Tungsten Carbide Tip With Lock 17 Cm- 6 3/4"

Needle Holder Mathieu-Kocher Tungsten Carbide Tip With Lock 17
Cm- 6 3/4"

[Read More](#)

SKU: BMN010C

Price:

Stock: instock

Categories: [Basic / General](#), [Micro Needle Holders](#)

Tags: [18.64.17](#)

Product Description

Needle Holder Mathieu-Kocher Tungsten Carbide Tip With Lock 17 Cm- 6 3/4"