

# Nelson Trocar 20 Cm-8" 4.0 Mm

Nelson Trocar 20 Cm-8" 4.0 Mm

[Read More](#)

**SKU:** UTR015

**Price:**

**Stock:** instock

**Categories:** [Trocar](#), [Urology](#)

**Tags:** [22.01.40](#)



22.01.22 2.2 mm

22.01.30 3.0 mm

22.01.40 4.0 mm

22.01.50 5.0 mm

22.01.55 5.5 mm

Nelson

## Product Description

Nelson Trocar 20 Cm-8" 4.0 Mm