



Osteotome 4mm Tip Working Length 30cms Fibre Handle Cvd

Osteotome 4mm Tip Working Length 30cms Fibre Handle Cvd

[Read More](#)

SKU: DMB527

Price:

Stock: instock

Categories: [Chisels](#), [Gouges](#), [Neuro & Spine](#), [Endoscopic Interlaminar](#), [Osteotome](#)

Tags: [DMB527](#), [FL556](#)

Product Description

Osteotome 4mm Tip Working Length 30cms Fibre Handle Cvd
