



Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 6mm

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 6mm

[Read More](#)

SKU: DMB293

Price:

Stock: instock

Categories: [Drills, Perforators, Trephines, Endoscopic Transilliac, Neuro & Spine, Perforators, Trephines](#)

Tags: [DMB293](#)

Product Description

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 6mm