



# Retractor Self Retaining Adson Hinged 26Cm 10 1/4" Sharp

Retractor Self Retaining Adson Hinged 26Cm 10 1/4" Sharp

[Read More](#)

**SKU:** NSR065

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Self Retaining Retractor](#)

**Tags:** [NSR065](#)

## Product Description

Retractor Self Retaining Adson Hinged 26Cm 10 1/4" Sharp

---