



# Retractor Self Retaining Beckmann-Adson 31 Cm 12 1/8" Sharp 25Mm

Retractor Self Retaining Beckmann-Adson 31 Cm 12 1/8" Sharp 25Mm

[Read More](#)

**SKU:** NSR100

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Self Retaining Retractor](#)

**Tags:** [31.24.31](#)

## Product Description

Retractor Self Retaining Beckmann-Adson 31 Cm 12 1/8" Sharp 25Mm