



Retractor Self Retaining Beckmann-Adson 31 Cm 12 1/8" Sharp 44Mm

Retractor Self Retaining Beckmann-Adson 31 Cm 12 1/8" Sharp 44Mm

[Read More](#)

SKU: NSR110

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Self Retaining Retractor](#)

Tags: [31.26.31](#)

Product Description

Retractor Self Retaining Beckmann-Adson 31 Cm 12 1/8" Sharp 44Mm