



Right Angle Ball Hook 5 Mm Tip Dia 3 Mm Bayonet 10"

Right Angle Ball Hook 5 Mm Tip Dia 3 Mm Bayonet 10"

[Read More](#)

SKU: MIS102

Price:

Stock: instock

Categories: [Hooks & Probes](#), [MIS Instruments](#), [Neuro & Spine](#)

Tags: [MIS102](#)

Product Description

Right Angle Ball Hook 5 Mm Tip Dia 3 Mm Bayonet 10"