



# Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 1.5mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 1.5mm

[Read More](#)

**SKU:** NSU085

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

**Tags:** [NSU085](#)

## Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 1.5mm