



# Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 1mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 1mm

[Read More](#)

**SKU:** NSU080

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

**Tags:** [NSU080](#)

## Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 1mm