



# Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 2.5mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working  
Length 10cm 2.5mm

[Read More](#)

**SKU:** NSU095

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

**Tags:** [NSU095](#)

## Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 10cm 2.5mm