



Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 13cm 2.5mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working
Length 13cm 2.5mm

[Read More](#)

SKU: NSU120

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU120](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 13cm 2.5mm