



Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 13cm 3mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 13cm 3mm

[Read More](#)

SKU: NSU125

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU125](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 13cm 3mm