



# Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 2.5mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 2.5mm

[Read More](#)

**SKU:** NSU145

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

**Tags:** [NSU145](#)

## Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 2.5mm