



Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 2mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 2mm

[Read More](#)

SKU: NSU140

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU140](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 2mm