



Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 10cm 1.5mm

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working
Length 10cm 1.5mm

[Read More](#)

SKU: NSU010

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU010](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 10cm 1.5mm