



Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 10cm 2.5mm

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working
Length 10cm 2.5mm

[Read More](#)

SKU: NSU020

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU020](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 10cm 2.5mm