



Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 13cm 1.5mm

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working
Length 13cm 1.5mm

[Read More](#)

SKU: NSU035

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU035](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 13cm 1.5mm