



Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 15cm 2mm

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working
Length 15cm 2mm

[Read More](#)

SKU: NSU065

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU065](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma Cvd Working Length 15cm 2mm