



Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 13cm 1.5mm

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 13cm 1.5mm

[Read More](#)

SKU: NSU185

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU185](#)

Product Description

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 13cm 1.5mm