



Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 15cm 2.5mm

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 15cm 2.5mm

[Read More](#)

SKU: NSU220

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU220](#)

Product Description

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 15cm 2.5mm