



Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 18cm 2.5mm

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working
Len 18cm 2.5mm

[Read More](#)

SKU: NSU245

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU245](#)

Product Description

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Len 18cm 2.5mm