



# Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Length 18cm 1mm

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working  
Length 18cm 1mm

[Read More](#)

**SKU:** NSU230

**Price:**

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

**Tags:** [NSU230](#)

## Product Description

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Length 18cm 1mm