



T Handle Cannula11.2mm Od 10.20mm Id Working Length 12cms Elevator Tip 10mm

T Handle Cannula11.2mm Od 10.20mm Id Working Length 12cms
Elevator Tip 10mm

[Read More](#)

SKU: DMB510

Price:

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB510](#), [FL525](#)

Product Description

T Handle Cannula11.2mm Od 10.20mm Id Working Length 12cms Elevator Tip 10mm