



# T Handle Sheath 17 Cms Chamfer 60 Deg L 15mm Id 6.5mm Od 7.5mm 6.3mm Chinese Scope

T Handle Sheath 17 Cms Chamfer 60 Deg L 15mm Id 6.5mm Od 7.5mm  
6.3mm Chinese Scope

[Read More](#)

**SKU:** DMB405C

**Price:**

**Stock:** instock

**Categories:** [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

**Tags:** [DMB405C](#), [FL224](#)

## Product Description

T Handle Sheath 17 Cms Chamfer 60 Deg L 15mm Id 6.5mm Od 7.5mm 6.3mm Chinese Scope