



T Handle Sheath 17 Cms Chamfer 60 Deg L 15mm Id 6.5mm Od 7.5mm 6.3mm Chinese Scope

T Handle Sheath 17 Cms Chamfer 60 Deg L 15mm Id 6.5mm Od 7.5mm
6.3mm Chinese Scope

[Read More](#)

SKU: DMB405C

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB405C](#), [FL224](#)

Product Description

T Handle Sheath 17 Cms Chamfer 60 Deg L 15mm Id 6.5mm Od 7.5mm 6.3mm Chinese Scope