



T Handle Sheath 17 Cms Elevator Tip 15mm 6.5mm Od 7.5mm 6.3mm Chinese Scope

T Handle Sheath 17 Cms Elevator Tip 15mm 6.5mm Od 7.5mm 6.3mm
Chinese Scope

[Read More](#)

SKU: DMB405D

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB405D](#), [FL226](#)

Product Description

T Handle Sheath 17 Cms Elevator Tip 15mm 6.5mm Od 7.5mm 6.3mm Chinese Scope