



T Handle Sheath 17Cms Chamfer 120 Deg L 12mm Id 6.5mm Od 7.5mm 6.3mm Chinese Scope

T Handle Sheath 17Cms Chamfer 120 Deg L 12mm Id 6.5mm Od
7.5mm 6.3mm Chinese Scope

[Read More](#)

SKU: DMB405B

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB405B](#), [FL187](#)

Product Description

T Handle Sheath 17Cms Chamfer 120 Deg L 12mm Id 6.5mm Od 7.5mm 6.3mm Chinese Scope