



T Handle Sheath 17Cms Chamfer 120 Deg L 12mmId 7mm Od 8mm 6.8mm EU Scope

T Handle Sheath 17Cms Chamfer 120 Deg L 12mmId 7mm Od 8mm
6.8mm EU Scope

[Read More](#)

SKU: DMB305B

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB305B](#), [FK507](#)

Product Description

T Handle Sheath 17Cms Chamfer 120 Deg L 12mmId 7mm Od 8mm 6.8mm EU Scope