



T Handle Sheath 17Cms Chamfer 60 Deg L 15mmId 7mm Od 8mm 6.8mm EU Scope

T Handle Sheath 17Cms Chamfer 60 Deg L 15mmId 7mm Od 8mm
6.8mm EU Scope

[Read More](#)

SKU: DMB305C

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB305C](#), [FK511](#)

Product Description

T Handle Sheath 17Cms Chamfer 60 Deg L 15mmId 7mm Od 8mm 6.8mm EU Scope