



T Handle Sheath 17Cms Elevator Tip 15mm Id 7mm Od 8mm (For 6.8mm EU Scope)

T Handle Sheath 17Cms Elevator Tip 15mm Id 7mm Od 8mm (For 6.8mm EU Scope)

[Read More](#)

SKU: DMB305D

Price:

Stock: instock

Categories: [Endoscopic Transforaminal](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB305D](#), [FK505](#)

Product Description

T Handle Sheath 17Cms Elevator Tip 15mm Id 7mm Od 8mm (For 6.8mm EU Scope)