



Tian Song Endoscope 30 Degree Angle 6.3mm Diameter & 3.5mm Work Channel HD Scope

Tian Song Endoscope 30 Degree Angle 6.3mm Diameter & 3.5mm Work Channel HD Scope

[Read More](#)

SKU: ESC100

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Endoscopic Transforaminal](#), [Endoscope](#)

Tags: [ESC100](#)

Product Description

Tian Song Endoscope 30 Degree Angle 6.3mm Diameter & 3.5mm Work Channel HD Scope