

Wire Twisters



70.01.15 15 cm - 6"
70.01.18 17 cm - 7"
70.01.20 20.3 cm - 8"

Wire Twister 15Cm-6" Tungsten Tip

Wire Twister 15Cm-6" Tungsten Tip

[Read More](#)

SKU: CWI005C

Price:

Stock: instock

Categories: [Cardio Thoracic](#), [Wire Instruments](#)

Tags: [70.01.15](#), [BM111](#)

Product Description

Wire Twister 15Cm-6" Tungsten Tip