

Wire Twisters



70.01.15 15 cm - 6"  
70.01.18 17 cm - 7"  
70.01.20 20.3 cm - 8"

# Wire Twister 20.3Cm-8" Tungsten Tip

Wire Twister 20.3Cm-8" Tungsten Tip

[Read More](#)

**SKU:** CWI015C

**Price:**

**Stock:** instock

**Categories:** [Cardio Thoracic](#), [Wire Instruments](#)

**Tags:** [70.01.20](#), [BM112](#)

## Product Description

Wire Twister 20.3Cm-8" Tungsten Tip